Cytotoxicity and cytokine expression induced by silorane and methacrylate-based composite resins

Daniele Lucca LONGO¹, Francisco Wanderley Garcia PAULA-SILVA², Lucia Helena FACCIOLI², Patrícia Maria GATÓN-HERNÁNDEZ^{1,3}, Alexandra Mussolino de QUEIROZ¹, Léa Assed Bezerra da SILVA¹

- 1- Universidade de São Paulo, Faculdade de Odontologia de Ribeirão Preto, Departamento de Clínica Infantil, Ribeirão Preto, SP, Brasil.
- 2- Universidade de São Paulo, Faculdade de Ciências Farmacêuticas de Ribeirão Preto, Departamento de Análises Clínicas, Ribeirão Preto, SP, Brasil.
- 3- Universitat de Barcelona, Facultat d'Odontologia, Department d'Odontopediatria, Barcelona, España.

Corresponding address: Léa Assed Bezerra da Silva - Faculdade de Odontologia de Ribeirão Preto, Universidade de São Paulo - Avenida do Café s/n -Monte Alegre - 14040-904 - Ribeirão Preto - SP - Brazil - Phone: +55-16-3315-3965 - Fax: +55-16-3315-4102. e-mail: lea@forp.usp.br

Submitted: October 5, 2015 - Modification: April 9, 2016 - Accepted: April 11, 2016

ABSTRACT

he successful use of composite resins in Dentistry depends on physicochemical properties, but also on the biological compatibility of resins, because of the close association between pulp and dentin. Objective: The aim of this study was to evaluate cytotoxicity and cytokine production induced by light-cured or non-light-cured methacrylate-based and silorane composite resins in RAW 264.7 macrophages. Material and Methods: Cells were stimulated with the extracts from light-cured or non-light-cured composite resins. After incubation for 24 h, cytotoxicity was assessed with the lactate dehydrogenase (LDH) and methyl thiazolyl tetrazolium (MTT) assays, and total protein was quantified using the Lowry method. TNF- α detection was examined with an enzyme-linked immunosorbent assay (ELISA) conducted with cell supernatants after cell stimulation for 6, 12, and 24 h. Data were analyzed using one-way analysis of variance (ANOVA) and Tukey's post hoc test (a=0.05). Results: Kalore™ and Filtek™ Silorane were cytotoxic with or without light curing (p<0.05) after 24 h of incubation. KaloreTM stimulated the early production of TNF- α in comparison with control (p<0.05), whereas FiltekTM Silorane did not affect TNF- α levels after 6 and 12 h (p>0.05). However, after 24 h Filtek™ Silorane inhibited the production of TNF- α (p<0.05). Conclusions: KaloreTM and FiltekTM Silorane were cytotoxic regardless of light curing. The extract obtained from Kalore™ after 15 days of incubation stimulated the production of TNF- α , unlike that obtained from FiltekTM Silorane.

Keywords: Dental materials. Composite resins. Cytotoxicity. Cytokines.

INTRODUCTION

The main components of composite resins include the organic matrix; filler particles; the bonding agent, which connects the filler to the organic matrix; the activator system, which initiates polymerization; pigments, which impart the compound with colors similar to those of teeth; and polymerization inhibitors, which increase the useful life and working time of the material^{1,16}. Composite resins are dental materials commonly used to restore the structural integrity and function of teeth affected by caries, erosion, fracture or attrition¹¹.

The effective application of resin materials

heavily depends on their physicochemical properties. However, clinical success ultimately depends on the biological compatibility of the resins, because of the close relation between pulp and dentin²⁷. The filler material does not seem to affect biological compatibility, which mostly depends on the organic components of composite resins⁵. These components include monomers, such as bisphenol A-glycidyl methacrylate (Bis-GMA), triethylene glycol dimethacrylate (TEGDMA), urethane dimethacrylate (UDMA), 2-hydroxyethyl methacrylate (HEMA), and bisphenol A-ethoxylated dimethacrylate (Bis-EMA)4,16. Monomers that are not polymerized are released from conventional methacrylate-based

composite resins and have been associated with genotoxicity^{5,8,9,22,25}, estrogenicity^{5,22,26}, changes in the immune system^{10,17}, hypersensitivity, cytotoxicity^{3,5,8,29}, and the production of reactive oxygen species^{3,21}.

The use of low shrinkage monomers and high molecular weight has improved significantly the composite resins, although mechanical and chemical issues still remain, specifically regarding the polymerization shrinkage effects^{15,18}. New methacrylate-based composite resins have been developed with changes in the composition, structure, and nature of polymerization, such as the resin Kalore™ (GC FUJI, Kasugai, Japan), which uses the high molecular weight Dupont monomer DX-511⁴. Conversely, the composite Filtek™ Silorane, brought to market by 3M ESPE (Seefeld, Germany), does not have a methacrylate monomer as its organic matrix, instead, it uses siloxane and oxirane molecules²⁸. To date, no studies have evaluated the cytotoxicity and the induction of pro-inflammatory molecules mediated by these composite resins in RAW 264.7 macrophages.

Therefore, this study assesses the cytotoxic effects and the production of cytokines induced by light-cured or non-light-cured methacrylate-based and silorane composite resins in RAW 264.7 macrophages.

MATERIAL AND METHODS

Extract preparation

This study was carried out according to the standards of the International Organization for Standardization (ISO) no. 10993-5:20097.

In the absence of artificial light, composite resins were removed from their tubes with a sterile spatula #1, placed on sterile paper, and weighed (PG 503-S, Mettler Toledo®; Toledo, Ohio, USA). Cell viability and cytokine production were evaluated after the exposure of the cells to the following experimental groups: KaloreTM and FiltekTM Silorane (light-cured and non-light-cured) (Table 1).

Light curing was performed using a lightemitting diode (LED) device (RADII-CAL SDI Limited; Bayswater, Victoria, Australia) with a light intensity of 1,200 mW/cm² and a wavelength of 430-480 nm. The polymerization process used for the composite resins followed the manufacturer's instructions, 20 seconds for each increment of composite resin (20 mg). For the groups of resins that were not light cured, the materials were kept in the dark until the moment of use.

To extract their components, composite resins were placed in 12-well plates and properly identified, and 3 mL of culture medium (DMEM) containing 500 μl of gentamicin (10 mg/mL; Gibco; Grand Island, NY, USA) and 5 mL of penicillin (100 $\mu g/mL$; Gibco; Grand Island, NY, USA) and streptomycin (100 $\mu g/mL$; Gibco; Grand Island, NY, USA) were added, so that the material was fully covered by this solution for 30 min. The solution was then removed with a serological pipette and each well of the plate was washed with 3 mL of PBS and filled with 3 mL of DMEM. Composite resins were weighed and the extracts 20-80 mg of resin per mL of DMEM were obtained after 15 d of incubation.

Cell culture

The RAW 264.7 murine macrophage cell line was obtained from the American Type Culture Collection (ATCC; Rockville, MD, USA). The cells were grown in DMEM medium supplemented with 10% fetal bovine serum and 1% gentamicin (DMEM-c). After the formation of a monolayer, cells were harvested with plastic cell scrapers and centrifuged at 1,500 rpm for 10 min at 10°C (Eppendorf; Hamburg, Hamburg, Germany). After centrifugation, supernatants were discarded and 10 mL of DMEM-c was added to each tube of cells. The total number of cells was counted and the viability was determined in a Neubauer chamber (BOECO Germany; Hamburg, Hamburg, Germany) using Trypan blue (Gibco; Grand Island, NY, USA). Cells were then plated in 96-well culture plates (Cell Wells - 25,820, Corning Glass Works; New York city, NY, USA) at a density of 1x10⁵ cells/ well and incubated overnight in DMEM-c in an incubator with a moist atmosphere of 5% CO₃ and 95% air at 37°C. Cells were then incubated with the extracts of composite resin for 24 h.

Table 1- Compositions of KaloreTM and Filtek™ Silorane resins, according to the information provided by the manufacturers

Composite Resin	Manufacturer/ Color	Particle	Organic Matrix	Inorganic Load (C.I.)	% C.I.
Kalore™	GC FUJI/A3	Nanohybrid	DX-511 UDMA dimethacrylatenonspecific comonomers	Glass Fluoroaminosilicate (silica and silicon dioxide)	82
Filtek™ Silorane	3M ESPE/A3	Microhybrid	Siloxane Oxirane	Quartz Yttrium fluoride	76

Lactate dehydrogenase (LDH) assay

The cytotoxicity of the resin extract in RAW 264.7 macrophages was evaluated through the level of LDH released in the cell supernatant following cell lysis using the CytoTox96® non-radioactive cytotoxicity assay (Promega Corporation; Madison, WI, USA). The absorbance was measured at 490 nm with a spectrophotometer (mQuanti, Bio-Tek Instruments, Inc.; Winooski, VT, USA). LDH levels were expressed as percentages of the LDH levels observed in control cultures.

MTT assay

Cell viability was evaluated using the 3-(4,5-dimethylthiazol-2-yl)-2,5-diphenyltetrazolium bromide (MTT) colorimetric assay (Sigma-Aldrich)²⁵. The cells were incubated with the extracts of composite resin for 24 h, then the cultures were incubated with 5% MTT in RPMI for 3 h. Subsequently, 50 mL of 20% sodium dodecyl sulfate (SDS) in 0.01 M HCl were added to each well and maintained at room temperature until the precipitate completely solubilized. Absorbance was measured at 570 nm with a spectrophotometer (mQuanti) and was directly proportional to cell viability. The cytotoxicity of the composite resins was expressed as percentages of the cytotoxicity observed in non-stimulated control cells.

Total protein quantification

Total protein quantification was performed using the Lowry method (Bio-rad DC Protein assay). Absorbance was measured at 750 nm with a spectrophotometer (mQuanti). Data are expressed as mg/mL obtained based on a standard curve using bovine serum albumin (BSA).

TNF- α detection

The concentration of TNF- α in culture supernatants was quantified by ELISA using specific purified and biotinylated antibodies and cytokine standards, according to the manufacturers' instructions (R & D Systems; Minneapolis, MN, USA). The optical densities were measured at 450 nm by a microplate reader. The cytokine concentrations were determined using a standard curve established with the appropriate recombinant cytokine, and are expressed in pg/mL.

Statistical analyses

Data represent the mean \pm SEM. Statistical variations were determined by one-way ANOVA and Tukey's test. Values of p<0.05 were considered to be significant.

RESULTS

Methacrylate-based and silorane composite resins are cytotoxic

Methacrylate-based and silorane composite resins induced cell lysis and LDH release into the cell supernatant regardless of whether resins were light cured or not (Figure 1A). Intracellular dehydrogenase activity, measured by MTT, was reduced by the methacrylate-based composite resin with or without light curing (p<0.05). A similar effect was observed when cells were exposed to the extract of non-light-cured silorane composite resin, but not with the light-cured resin (Figure 1B). Total protein quantification indicated that methacrylate-based and silorane composite resins inhibited cell activity (p<0.05, Figure 1C).




Figure 1- Effects of different composite resin extracts (Kalore[™] and Filtek[™] Silorane, light-cured and non-light-cured) on murine RAW 264.7 macrophage cultures, as shown by the LDH assay (A), MTT assay (B), and total protein quantification (C), after 24-hour incubation. Statistically significant differences (p<0.05) in relation to the control group (macrophages incubated in culture medium only) are indicated by *




Figure 2- Effects of different composite resin extracts (KaloreTM and FiltekTM Silorane, light-cured and non-light-cured) on the release of TNF- α by cultured murine RAW 264.7 macrophages, as assessed by ELISA after incubation periods of 6 (A, B), 12 (C, D), and 24 (E, F) h. Statistically significant differences (p<0.05) in relation to the control group (macrophages incubated in culture medium only) are indicated by *

Macrophage TNF- α production is stimulated by methacrylate-based composite resin and inhibited by silorane

The extract of the methacrylate-based composite resin stimulated the early production of TNF- α in comparison with the control treatment (p<0.05, Figure 2A and 2C). The silorane composite resin extract did not elicit statistically significant differences in the production of TNF- α relative to the control after 6 and 12 h of incubation (p>0.05, Figure 2B and 2D). Interestingly, after 24 h of incubation, the methacrylate-based composite resin sustained TNF- α release whereas the silorane resin extract inhibited the production of TNF- α compared with the control (p<0.05, Figure 2E and 2F).

DISCUSSION

Methacrylate-based and silorane composite resins were cytotoxic with or without light curing as indicated by LDH release and total protein quantification. Extracts of methacrylate-based composite resins stimulated the production of TNF- α after 6 h of incubation. On the other hand, silorane resin extract inhibited the production of TNF- α after 24 h of incubation.

Non-light-curing or incomplete light-curing of methacrylate-based composites results in the release of resin matrix components, called residual monomers, i.e., unpolymerized monomers^{2,21}. The release of these non-light-cured resin materials has been associated with various adverse effects^{3,5,8-10,17,21,22,25,26,29}. Regarding cytotoxic effects specifically, it is known that isolated monomers cause several biological effects on cells, such as damage to the cell membrane, inhibition of metabolic enzyme activity, cell-cycle delay and interruption, gene mutation, DNA breakage and apoptosis after reduction of GSH (glutathione) via oxidative stress²². Several previous studies have investigated the toxicity of monomers in isolation in various cell types^{3,5,8,22,29}. However, a recent study using RAW 264.7 macrophages showed that co-exposure to TEGDMA and both types of filler particles, Nanosilica and Quartz, resulted in an additive attenuation of the LPS-induced IL-1B release, but cellular viability and TNF- α release were not significantly affected¹³. On the other hand, the assessment of these isolated monomers does not faithfully reproduce the procedure for the routine use of these materials in the clinic, where they are applied in the polymerized form, i.e., in the form of composite resin.

Our findings revealed that the methacrylatebased composite was cytotoxic *in vitro* in RAW 264.7 mouse macrophages. Similarly, another study showed the biological effects of metacrylate-based composite resins, including possible changes to the DNA of skin cells in class V restorations²³, and the authors concluded that the resin may cause cellular damage in fibroblasts²⁴. However, the silorane composite resin did not affect intracellular mitochondrial dehydrogenase activity, except when it was not light-cured. The difference between the observed cytotoxicity in this study, and the satisfactory compatibility obtained in other studies by our group^{19,20} may result from complete or incomplete curing of the composite resin, respectively. Although we observed a cytotoxic effect of the silorane composite resin, this was not reflected on the MTT assay. Likewise, another study suggested the non-toxic nature of the silorane composite resin in human fibroblasts treated with Filtek P90, which showed only an insignificant decrease in cell proliferation in 24 h and 48 h¹².

The macrophage lineage, like many cells of the immune system, plays a role in innate immunity-related functions and in the production of inflammatory mediators, such as TNF- α^{14} . The metacrylate-based composite resin sustained TNF- α production for up to 24 h of incubation (80 mg/mL), although lower concentrations of the extracts (20-40 mg/mL) stimulated only the early production of TNF- α (6 h). This contrast can be explained by cell death, which decreases the production of TNF- α .

Finally, the silorane composite resin did not induce TNF- α after 6 and 12 h, but inhibited it after 24 h. The decrease in TNF- α levels induced by unpolimeryzed silorane composite resin can also be associated with cytotoxicity. However, that reduction also happened when the resin was polymerized, and in parallel with unchanged cell viability. Thus, the silorane composite resin induced an anti-inflammatory response. Another study revealed that the silorane composite resin had a negative effect on TNF- α levels⁶.

We conclude that KaloreTM and FiltekTM Silorane were cytotoxic regardless of light curing. Interestingly, KaloreTM stimulated the production of TNF- α , unlike FiltekTM Silorane.

ACKNOWLEDGEMENTS

The authors deny any conflicts of interest related to this study. This study was partially supported by grants from the São Paulo Research Foundation (FAPESP 2010/17611-4 to FWGPS) and CNPq (135819/2011-3 fellowship to DLL).

REFERENCES

- 1- Asmussen E, Peutzfeldt A. Influence of UEDMA BisGMA and TEGDMA on selected mechanical properties of experimental resin composites. Dent Mater. 1998;14(1):51-6.
- 2- Brackett MG, Brackett WW, Browning WD, Rueggeberg FA. The effect of light curing source on the residual yellowing of resin composites. Oper Dent. 2007;32(5):443-50.

- 3- Chang MC, Chen LI, Chan CP, Lee JJ, Wang TM, Yang TT, et al. The role of reactive oxygen species and hemeoxygenase-1 expression in the cytotoxicity, cell cycle alteration and apoptosis of dental pulp cells induced by BisGMA. Biomaterials. 2010;31(32):8164-71.
- 4- Ferracane JL. Resin composite state of the art. Dent Mater. 2011;27(1):29-38.
- 5- Goldberg M. *In vitro* and *in vivo* studies on the toxicity of dental resin components: a review. Clin Oral Investig. 2008;12(1):1-8. 6- Ilday NO, Celik N, Dilsiz A, Alp HH, Aydin T, Seven N, et al. The effects of silorane composites on levels of cytokines and
- periodontal parameters. Contemp Clin Dent. 2013;4(4):437-42.
 7- International Organization for Standardization. ISO 10993-
- 5: Biological evaluation of medical devices Tests for *in vitro* cytotoxicity. Geneva: ISO; 2009.
- 8- Kleinsasser NH, Schmid K, Sassen AW, Harréus UA, Staudenmaier R, Folwaczny M, et al. Cytotoxic and genotoxic effects of resin monomers in human salivary gland tissue and lymphocytes as assessed by the single cell microgel electrophoresis (Comet) assay. Biomaterials. 2006;27(9):1762-70.
- 9- Kleinsasser NH, Wallner BC, Harréus UA, Kleinjung T, Folwaczny M, Hickel R, et al. Genotoxicity and cytotoxicity of dental materials in human lymphocytes as assessed by the single cell microgel electrophoresis (comet) assay. J Dent. 2004;32(3):229-34.
- 10- Kostoryz EL, Tong PY, Strautman AF, Glaros AG, Eick JD, Yourtee DM. Effects of dental resins on TNF-alpha-induced ICAM-1 expression in endothelial cells. J Dent Res. 2009;80(9):1789-92.
- 11- Krifka S, Spagnuolo G, Schmalz G, Schweikl H. A review of adaptive mechanisms in cell responses towards oxidative stress caused by dental resin monomers. Biomaterials. 2013;34(19):4555-63.
- 12- Madhyastha PS, Naik DG, Kotian R, Padma D, Srikant N, Bhat KM. Evaluation of cytotoxicity of silorane and methacrylate based dental composites using human gingival fibroblasts. J Clin Diagn Res. 2015;9(1):5-8.
- 13- Mathisen GH, Ansteinsson V, Samuelsen JT, Becher R, Dahl JE, Bølling AK. TEGDMA and filler particles from dental composites additively attenuate LPS-induced cytokine release from the macrophage cell line RAW 264.7. Clin Oral Investig. 2015;19(1):61-9.
- 14- Metzger Z. Macrophages in periapical lesions. Endod Dent Traumatol. 2000;16(1):1-8.
- 15- Naoum SJ, Ellakwa A, Morgan L, White K, Martin FE, Lee IB. Polymerization profile analysis of resin composite dental restorative materials in real time. J Dent. 2012;40(1):64-70.

- 16- Peutzfeldt A. Resin composites in dentistry: the monomer systems. Eur J Oral Sci. 1997;105(2):97-116.
- 17- Rakich DR, Wataha JC, Lefebvre CA, Weller RN. Effect of dentin bonding agents on the secretion of inflammatory mediators from macrophages. J Endod. 1999;25(2):114-7.
- 18- Rodrigues AS Jr, Scherrer SS, Ferracane JL, Della Bona A. Microstructural characterization and fracture behavior of a microhybrid and a nanofil composite. Dent Mater. 2008;24:1281-8. 19- Ruiz-de-Castañeda E, Gatón-Hernández P, Rodriguez EG, Silva RA, Nelson-Filho P, Silva LA. Pulpal and periapical response after restoration of deep cavities in dogs' teeth with Filtek Silorane and Filtek Supreme XT systems. Oper Dent. 2013;38(1):73-81. 20- Ruiz-de-Castañeda E, Silva LA, Gaton-Hernández P, Consolaro A, Rodriguez EG, Silva RA, et al. Filtek™ Silorane and Filtek™ Supreme XT resins: tissue reaction after subcutaneous implantation in isogenic mice. Braz Dent J. 2011;22(2):105-10.
- dental monomers. Adv Dent Res. 2011;23(3):302-6. 22- Schweikl H, Hiller KA, Bolay C, Kreissl M, Kreismann W, Nusser A, et al. Cytotoxic and mutagenic effects of dental composite materials. Biomaterials. 2005;26:1713-9.

21- Schmalz G, Krifka S, Schweikl H. Toll-like receptors, LPS, and

- 23- Tadin A, Galic N, Mladinic M, Marovic D, Kovacic I, Zeljezic D. Genotoxicity in gingival cells of patients undergoing tooth restoration with two different dental composite materials. Clin Oral Investig. 2014;18(1):87-96.
- 24- Tadin A, Marovic D, Galic N, Kovacic I, Zeljezic D. Composite-induced toxicity in human gingival and pulp fibroblast cells. Acta Odontol Scand. 2014;72(4):304-11.
- 25- Urcan E, Scherthan H, Styllou M, Haertel U, Hickel R, Reichl FX. Induction of DNA double-strand breaks in primary gingival fibroblasts by exposure to dental resin composites. Biomaterials. 2010;31(8):2010-4.
- 26- Wada H, Turami H, Imazato S, Narimatsu M, Ebisu S. *In vitro* estrogenicity of resin composites. J Dent Res. 2004;83:222-6.
- 27- Wataha JC. Predicting clinical biological responses to dental materials. Dent Mater. 2012;28(1):23-40.
- 28- Weinmann W, Thalacker C, Guggenberger R. Siloranes in dental composites. Dent Mater. 2005;21:68-74.
- 29- Wisniewska-Jarosinska M, Poplawski T, Chojnacki CJ, Pawlowska E, Krupa R, Szczepanska J, et al. Independent and combined cytotoxicity and genotoxicity of triethylene glycol dimethacrylate and urethane dimethacrylate. Mol Biol Rep. 2011;38(7):4603-11.